
Wir verstehen Hören.

bruckhoff

PRODUCT CATALOGUE

BRUCKHOFF HEARING SYSTEMS PORTFOLIO.



- 2 <u>BRUCKHOFF</u>
 - • •
 - • •
 - . . .
 - • •
 - • •
 - • •
 - • •
 - • •
 - • •
 - • •









THE SPECIALIST FOR BONE CONDUCTION.

Bruckhoff develops and produces technically and optically innovative Hearing Systems and Medicine Products. To the best known own Products belong the Hearing Glasses La Belle as well as Products of Bone Conduction Hearing Systems.



JUNIOR BONE CONDUCTION. CHILDRENS HEARING SYSTEM.

HEADBANDS WITH BONE CONDUCTION HEARING SYSTEM IN DIFFERENT COLOURS.



SO LIGHT, SO GOOD.

Carefree hearing, carefree playing: The combination of headband and hearing module offers many advantages to your child. Thus, the junior children's hearing systems is

- · hardly noticeable for your child,
- · extremely comfortable to wear,
- · automatically positioned correctly on the head,
- · almost invisible to others.

A children's hearing system of high wearing comfort: pleasant, comfortable and discreet. The junior hearing system is made in Germany. A quality product that will accord your child long-lasting pleasure.

SMALL BUT WOW!

- · barely visible from the outside
- adjustable headband
- simple and quick adjustment
- no surgery necessary
- good sound quality
- noise reduction

FASHIONABLE HEADBANDS

The customers receive a great selection of beautiful and fashionable headbands for the child. These are already prefabricated for easy assembly of the junior hearing system.

Even personal headbands can be modified for deployment of the junior hearing system.



- wind noise reduction
- two levels of feedback canceling
- phonak Roger system usable
- volume control and program selection can be deactivated
- fitting with bone conduction fitting software
- compatible with NOAH and NOAH Link from Vers. 4.3

LA BELLE BONE CONDUCTION. SPECTACLE HEARING SYSTEM.

UNIQUE COMBINATION OF GLASSES AND HEARING MODULE.



INGENIOUS: THE GLASSES BECOME THE HEARING SYSTEM.

The patented hearing module is technically and optically integrated into the glasses so that the added value of the glasses is virtually invisible. Due to the familiar position of the glasses on the head, the new hearing system is

- · automatically ideally positioned on the head,
- not noticeable to the wearer,
- hardly visible to others.

A hearing system of high wearing comfort: pleasant, comfortable and discreet. The la belle hearing system is made in Germany. A quality product that will guarantee the wearer long-lasting pleasure.




Practical. Whether you wear work, reading or sunglasses, you only need one hearing module. Click on the glasses and you're done!

TECHNOLOGY AT ITS BEST.

- · realistic sound
- · automatic adjustment of the listening intelligibility
- good hearing in noisy surroundings
- noise reduction
- wind noise reduction
- two levels of feedback canceling
- fitting with bone conduction fitting software

- compatible with NOAH and NOAH Link from Vers. 4.3
- end pieces cold bendable
- · this makes the position of the receiver flexible
- chamber lid in gold, silver or anthracite possible
- battery compartment cover size 13
- change to different glasses possible
- · almost all glasses can be used for the click system





HEADSET BONE CONDUCTION. HEADSET HEARING SYSTEM.

INNOVATIVE HEADSET WHICH ADAPTS TO ANY HEAD SHAPE.



PUT ON AND DONE

The patented hearing module is mounted on an adjustable metal bracket and resembles a headset or headphones. Due to the predefined position of the hearing modules the new hearing system is

- · ideally positioned on the head in a jiffy,
- light and thus hardly noticeable.

The headset hearing system is made in Germany. A quality product that will guarantee the wearer long-lasting pleasure.

TECHNOLOGY AT ITS BEST.

- easy to use
- realistic sound
- · automatic adjustment of the listening intelligibility
- perfect hearing in noisy surroundings
- noise reduction
- wind noise reduction
- two levels of feedback canceling

FROM SMALL TO LARGE.

The light and subtle metal bracket can be quickly and easily adjusted to different head sizes via the two adjusters. In addition, the bend of medical spring steel can be changed. Thus, the hearing system can be placed in the perfect position.



- fitting with bone conduction fitting software
- compatible with NOAH and NOAH Link from Vers. 4.3
- end pieces cold bendable
- this makes the position of the receiver flexible
- chamber lid in gold, silver or anthracite possible
- battery compartment cover size 13
- adjustable headband

BC-FIT SOFTWARE. MANAGE EASILY.

FITTING SOFTWARE FOR ALL BRUCKHOFF PRODUCTS.


FITTINGS SOFTWARE AT A GLANCE.

The BC-Fit fitting software provides a tool for trouble-free management of all your day-to-day tasks. BC-Fit is suitable for:

- customer management
- device selection, programming and adjustment
- providing advice (data logging / recommending adjustments)
- an overview of information about the hearing system

The optimum position of the earpiece is determined with the help of a tone generator. A 'first fit' is also possible at the touch of a button.

Download our BC-Fit fitting software now from www.bruckhoff. com in the "Downloads" section.



One-Click. Configure and adapt your customers' hearing systems at the touch of a button using pre-defined profiles.



Everything at a glance. The 'Data history' section provides an insight into the customer's habits, e.g. how often the volume has been changed using the VC adjuster or how often programme 2 has been used etc. The display shows an extract of the data history.



Practical. Clear display allows the right and left side of the hearing system to be easily and intuitively adjusted.

MATRIX MODULES LA BELLE BC BONE CONDUCTION

	LA BELLE BC EYEGLASSES		LA BELLE JUNIOR BC Headband	LA BELLE HEADSET BC Headset
FEATURES	LA BELLE BC D50	LA BELLE BC D70	LA BELLE JUNIOR BC D50	LA BELLE HEADSET BC D50
WDRC	16 channels	16 channels	16 channels	16 channels
Feedback canceler	4. generation	4. generation	4. generation	4. generation
Layered noise reduction™	•	•	•	•
Wind noise suppression	•	•	•	•
МРО	•	•	•	•
Volume control	•	•	•	•
Channels	16	16	16	16
In situ tone generator	•	•	•	•
Programmes	4	4	4	4
Push button	programmes	programmes	programmes	programmes
Program switch tones	•	•	•	•
Low battery tone	•	•	•	•
Telecoil	•	•	•	•
Microphones	1	2	1	1
Adaptive directional		•		
Audio input	•	•	•	•
Phonak Roger X receiver, Option			•	
Battery size 13	•	•	•	•
Fitting software special for BC	BC-Fit	BC-Fit	BC-Fit	BC-Fit
NOAH Vers. 4.3 and above	•	•	•	•
NOAH Link	•	•	•	•